

MRAS

Reactivity: Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 24kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MRAS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

M-RAS; RRAS3; R-RAS3;

Catalog #: A8315

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 22808

Isotype: IgG

Swiss Prot: O14807

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the Ras family of small GTPases. These membrane-associated proteins function as signal transducers in multiple processes including cell growth and differentiation, and dysregulation of Ras signaling has been associated with many types of cancer. The encoded protein may play a role in the tumor necrosis factor-alpha and MAP kinase signaling pathways. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

To place an order, please [Click HERE](#).