

DCTN6

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 21kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human DCTN6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

WS3; p27; WS-3;

Catalog #: A8316

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10671

Isotype: IgG

Swiss Prot: O00399

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene contains an RGD (Arg-Gly-Asp) motif in the N-terminal region, which confers adhesive properties to macromolecular proteins like fibronectin. It shares a high degree of sequence similarity with the mouse homolog, which has been suggested to play a role in mitochondrial biogenesis. The exact biological function of this gene is not known.

To place an order, please [Click HERE](#).