

SLC9A9

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 73kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SLC9A9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NHE9; AUTS16;

Catalog #: A8323

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 285195

Isotype: IgG

Swiss Prot: Q8IVB4

Purity: Affinity purification

For research use only.

Background:

This gene encodes a sodium/proton exchanger that is a member of the solute carrier 9 protein family. The encoded protein localizes to the late recycling endosomes and may play an important role in maintaining cation homeostasis. Mutations in this gene are associated with autism susceptibility 16 and attention-deficit/hyperactivity disorder.

To place an order, please [Click HERE](#).