

MID2

Reactivity:Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:83kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human MID2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FXY2; RNF60; TRIM1; MRX101;

Catalog #:A8325

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:11043

Isotype:IgG

Swiss Prot:Q9UJV3

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein localizes to microtubular structures in the cytoplasm. Alternate splicing of this gene results in two transcript variants encoding different isoforms.

To place an order, please [Click HERE](#).