

RNF11

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 17kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RNF11

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CGI-123; SID1669;

Catalog #: A8328

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 26994

Isotype: IgG

Swiss Prot: Q9Y3C5

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene contains a RING-H2 finger motif, which is known to be important for protein-protein interactions. The expression of this gene has been shown to be induced by mutant RET proteins (MEN2A/MEN2B). The germline mutations in RET gene are known to be responsible for the development of multiple endocrine neoplasia (MEN).

To place an order, please [Click HERE](#).