

UBE2D4

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 17kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human UBE2D4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HBUCE1;

To place an order, please [Click HERE](#).

Catalog #: A8332

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51619

Isotype: IgG

Swiss Prot: Q9Y2X8

Purity: Affinity purification

For research use only.