

SMAD7

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 46kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SMAD7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CRCS3; MADH7; MADH8;

Catalog #: A8342

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4092

Isotype: IgG

Swiss Prot: O15105

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a nuclear protein that binds the E3 ubiquitin ligase SMURF2. Upon binding, this complex translocates to the cytoplasm, where it interacts with TGF-beta receptor type-1 (TGFBR1), leading to the degradation of both the encoded protein and TGFBR1. Expression of this gene is induced by TGFBR1. Variations in this gene are a cause of susceptibility to colorectal cancer type 3 (CRCS3). Several transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).