

MYBPC1

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 128kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MYBPC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LCCS4; MYBPCC; MYBPCS;

Catalog #: A8348

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4604

Isotype: IgG

Swiss Prot: Q00872

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the myosin-binding protein C family. Myosin-binding protein C family members are myosin-associated proteins found in the cross-bridge-bearing zone (C region) of A bands in striated muscle. The encoded protein is the slow skeletal muscle isoform of myosin-binding protein C and plays an important role in muscle contraction by recruiting muscle-type creatine kinase to myosin filaments. Mutations in this gene are associated with distal arthrogryposis type I. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

To place an order, please [Click HERE](#).