

NME4

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 21kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NME4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NDPK-D; NM23H4; nm23-H4;

Catalog #: A8350

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4833

Isotype: IgG

Swiss Prot: O00746

Purity: Affinity purification

For research use only.

Background:

The nucleoside diphosphate (NDP) kinases (EC 2.7.4.6) are ubiquitous enzymes that catalyze transfer of gamma-phosphates, via a phosphohistidine intermediate, between nucleoside and dioxynucleoside tri- and diphosphates. The enzymes are products of the nm23 gene family, which includes NME4 (Milon et al., 1997)

To place an order, please [Click HERE](#).