

RNF5

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 20kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RNF5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RMA1; RING5;

Catalog #: A8351

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6048

Isotype: IgG

Swiss Prot: Q99942

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene contains a RING finger, which is a motif known to be involved in protein-protein interactions. This protein is a membrane-bound ubiquitin ligase. It can regulate cell motility by targeting paxillin ubiquitination and altering the distribution and localization of paxillin in cytoplasm and cell focal adhesions.

To place an order, please [Click HERE](#).