

RPN2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 69kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RPN2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SWP1; RPNII; RIBIIR; RPN-II;

Catalog #: A8352

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6185

Isotype: IgG

Swiss Prot: P04844

Purity: Affinity purification

For research use only.

Background:

This gene encodes a type I integral membrane protein found only in the rough endoplasmic reticulum. The encoded protein is part of an N-oligosaccharyl transferase complex that links high mannose oligosaccharides to asparagine residues found in the Asn-X-Ser/Thr consensus motif of nascent polypeptide chains. This protein is similar in sequence to the yeast oligosaccharyl transferase subunit SWP1. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).