

HLTF

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 114kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human HLTF

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ZBU1; HLTF1; RNF80; HIP116; SNF2L3; HIP116A; SMARCA3;

Catalog #: A8353

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6596

Isotype: IgG

Swiss Prot: Q14527

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the SWI/SNF family. Members of this family have helicase and ATPase activities and are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein contains a RING finger DNA binding motif. Two transcript variants encoding the same protein have been found for this gene. However, use of an alternative translation start site produces an isoform that is truncated at the N-terminus compared to the full-length protein.

To place an order, please [Click HERE](#).