

PPP6R2

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 105kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PPP6R2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PP6R2; SAPS2; SAP190; KIAA0685;

Catalog #: A8359

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9701

Isotype: IgG

Swiss Prot: O75170

Purity: Affinity purification

For research use only.

Background:

Protein phosphatase regulatory subunits, such as SAPS2, modulate the activity of protein phosphatase catalytic subunits by restricting substrate specificity, recruiting substrates, and determining the intracellular localization of the holoenzyme. SAPS2 is a regulatory subunit for the protein phosphatase-6 catalytic subunit (PPP6C; MIM 612725) (Stefansson and Brautigan, 2006

To place an order, please [Click HERE](#).