

RNF13

Reactivity:Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:43kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RNF13

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RZF;

Catalog #:A8363

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:11342

Isotype:IgG

Swiss Prot:O43567

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene contains a RING zinc finger, a motif known to be involved in protein-protein interactions. The specific function of this gene has not yet been determined.

Alternatively spliced transcript variants that encode the same protein have been reported. A pseudogene, which is also located on chromosome 3, has been defined for this gene.

To place an order, please [Click HERE](#).