

MSRB2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 20kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MSRB2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CBS1; MSRB; PILB; CBS-1; CGI-131;

To place an order, please [Click HERE](#).

Catalog #: A8364

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 22921

Isotype: IgG

Swiss Prot: Q9Y3D2

Purity: Affinity purification

For research use only.