

DHTKD1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 103kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human DHTKD1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CMT2Q; AMOXAD;

Catalog #: A8369

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55526

Isotype: IgG

Swiss Prot: Q96HY7

Purity: Affinity purification

For research use only.

Background:

This gene encodes a component of a mitochondrial 2-oxoglutarate-dehydrogenase-complex-like protein involved in the degradation pathways of several amino acids, including lysine. Mutations in this gene are associated with 2-aminoadipic 2-oxoadipic aciduria and Charcot-Marie-Tooth Disease Type 2Q.

To place an order, please [Click HERE](#).