

METTL3

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 64kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human METTL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

M6A; IME4; Spo8; MT-A70;

Catalog #: A8370

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 56339

Isotype: IgG

Swiss Prot: Q86U44

Purity: Affinity purification

For research use only.

Background:

This gene encodes the 70 kDa subunit of MT-A which is part of N6-adenosine-methyltransferase.

This enzyme is involved in the posttranscriptional methylation of internal adenosine residues in eukaryotic mRNAs, forming N6-methyladenosine.

To place an order, please [Click HERE](#).