

BTLA

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:4000

Calculated MW: 33kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human BTLA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BTLA1; CD272;

Catalog #: A8377

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 151888

Isotype: IgG

Swiss Prot: Q7Z6A9

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the immunoglobulin superfamily. The encoded protein contains a single immunoglobulin (Ig) domain and is a receptor that relays inhibitory signals to suppress the immune response. Alternative splicing results in multiple transcript variants. Polymorphisms in this gene have been associated with an increased risk of rheumatoid arthritis.

To place an order, please [Click HERE](#).