

IL4I1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 63kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IL4I1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LAO; FIG1; LAAO;

Catalog #: A8378

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 259307

Isotype: IgG

Swiss Prot: Q96RQ9

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein with limited similarity to L-amino acid oxidase which contains the conserved amino acids thought to be involved in catalysis and binding of flavin adenine dinucleotide (FAD) cofactor. The expression of this gene can be induced by interleukin 4 in B cells, however, expression of transcripts containing the first two exons of the upstream gene is found in other cell types. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).