

TMEM189

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 31kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human TMEM189

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KUA;

Catalog #: A8380

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 387521

Isotype: IgG

Swiss Prot: A5PLL7

Purity: Affinity purification

For research use only.

Background:

Co-transcription of this gene and the neighboring downstream gene (ubiquitin-conjugating enzyme E2 variant 1) generates a rare read-through transcript, which encodes a fusion protein comprised of sequence sharing identity with each individual gene product. The protein encoded by this individual gene lacks a UEV1 domain but includes three transmembrane regions. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).