

ATAD2

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 159kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human ATAD2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ANCCA; ATAD2; DKFZp667N1320; PRO2000; L16

Catalog #: A8382

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 29028

Isotype: IgG

Swiss Prot: Q6PL18

Purity: Affinity purification

For research use only.

Background:

ATAD2 is one of them encoding a conserved factor harbouring an AAA type ATPase domain and a bromodomain. ATAD2 controls chromatin dynamics, genome transcriptional activities and apoptotic cell response. This protein has 2 isoforms produced by alternative splicing.

To place an order, please [Click HERE](#).