

GLI1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 118kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GLI1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GLI;

Catalog #: A8387

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2735

Isotype: IgG

Swiss Prot: P08151

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the Kruppel family of zinc finger proteins. The encoded transcription factor is activated by the sonic hedgehog signal transduction cascade and regulates stem cell proliferation. The activity and nuclear localization of this protein is negatively regulated by p53 in an inhibitory loop. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).