

RSPO3

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 31kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human RSPO3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PWTSR; THSD2; CRISTIN1;

Catalog #: A8389

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84870

Isotype: IgG

Swiss Prot: Q9BXY4

Purity: Affinity purification

For research use only.

Background:

This gene belongs to the R-spondin family. The encoded protein plays a role in the regulation of Wnt (wingless-type MMTV integration site family)/beta-catenin and Wnt/planar cell polarity (PCP) signaling pathways, which are involved in development, cell growth and disease pathogenesis. Genome-wide association studies suggest a correlation of this gene with bone mineral density and risk of fracture. This gene may be involved in tumor development.

To place an order, please [Click HERE](#).