

YTHDF3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 64kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human YTHDF3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A8395

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 253943

Isotype: IgG

Swiss Prot: Q7Z739

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the YTH (YT521-B homology) domain protein family. The YTH domain is common in eukaryotes, is often found in the middle of the protein sequence, and may function in binding to RNA. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).