

## AGXT

---

**Reactivity:** Human

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 43kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human AGXT

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

AGT; PH1; SPT; AGT1; SPAT; TLH6; AGXT1;

**Catalog #:** A8397

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 189

**Isotype:** IgG

**Swiss Prot:** P21549

**Purity:** Affinity purification

For research use only.

**Background:**

This gene is expressed only in the liver and the encoded protein is localized mostly in the peroxisomes, where it is involved in glyoxylate detoxification. Mutations in this gene, some of which alter subcellular targeting, have been associated with type I primary hyperoxaluria.

*To place an order, please [Click HERE](#).*