

ALAD

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 36kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ALAD

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PBGS; ALADH;

Catalog #: A8398

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 210

Isotype: IgG

Swiss Prot: P13716

Purity: Affinity purification

For research use only.

Background:

The ALAD enzyme is composed of 8 identical subunits and catalyzes the condensation of 2 molecules of delta-aminolevulinate to form uroporphobilinogen (a precursor of heme, cytochromes and other hemoproteins). ALAD catalyzes the second step in the porphyrin and heme biosynthetic pathway; zinc is essential for enzymatic activity. ALAD enzymatic activity is inhibited by lead and a defect in the ALAD structural gene can cause increased sensitivity to lead poisoning and acute hepatic porphyria.

To place an order, please [Click HERE](#).