

CLTB

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 25kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CLTB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LCB;

Catalog #: A8404

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1212

Isotype: IgG

Swiss Prot: P09497

Purity: Affinity purification

For research use only.

Background:

Clathrin is a large, soluble protein composed of heavy and light chains. It functions as the main structural component of the lattice-type cytoplasmic face of coated pits and vesicles which entrap specific macromolecules during receptor-mediated endocytosis. This gene encodes one of two clathrin light chain proteins which are believed to function as regulatory elements. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).