

CR2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 113kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CR2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CR; C3DR; CD21; CVID7; SLEB9;

Catalog #: A8407

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1380

Isotype: IgG

Swiss Prot: P20023

Purity: Affinity purification

For research use only.

Background:

This gene encodes a membrane protein, which functions as a receptor for Epstein-Barr virus (EBV) binding on B and T lymphocytes. Genetic variations in this gene are associated with susceptibility to systemic lupus erythematosus type 9 (SLEB9). Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).