

POLL

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:63kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human POLL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BETAN; POLKAPPA;

Catalog #:A8408

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:27343

Isotype:IgG

Swiss Prot:Q9UGP5

Purity:Affinity purification

For research use only.

Background:

This gene encodes a DNA polymerase. DNA polymerases catalyze DNA-template-directed extension of the 3'-end of a DNA strand. This particular polymerase, which is a member of the X family of DNA polymerases, likely plays a role in non-homologous end joining and other DNA repair processes. Alternatively spliced transcript variants have been described.

To place an order, please [Click HERE](#).