

ESRRB

Reactivity:Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:56kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ESRRB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ERR2; ERRb; ESRL2; NR3B2; DFNB35;

Catalog #:A8416

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2103

Isotype:IgG

Swiss Prot:O95718

Purity:Affinity purification

For research use only.

Background:

This gene encodes a protein with similarity to the estrogen receptor. Its function is unknown; however, a similar protein in mouse plays an essential role in placental development.

To place an order, please [Click HERE](#).