

GHRHR

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 47kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GHRHR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GRFR; GHRFR; IGHD1B;

Catalog #: A8421

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2692

Isotype: IgG

Swiss Prot: Q02643

Purity: Affinity purification

For research use only.

Background:

This gene encodes a receptor for growth hormone-releasing hormone. Binding of this hormone to the receptor leads to synthesis and release of growth hormone. Mutations in this gene have been associated with isolated growth hormone deficiency (IGHD), also known as Dwarfism of Sindh, a disorder characterized by short stature.

To place an order, please [Click HERE](#).