

GNB2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 37kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GNB2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A8422

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2783

Isotype: IgG

Swiss Prot: P62879

Purity: Affinity purification

For research use only.

Background:

Heterotrimeric guanine nucleotide-binding proteins (G proteins), which integrate signals between receptors and effector proteins, are composed of an alpha, a beta, and a gamma subunit. These subunits are encoded by families of related genes. This gene encodes a beta subunit. Beta subunits are important regulators of alpha subunits, as well as of certain signal transduction receptors and effectors. This gene contains a trinucleotide (CCG) repeat length polymorphism in its 5' UTR.

To place an order, please [Click HERE](#).