

GNRH2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 13kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GNRH2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GnRH-II; LH-RHII;

Catalog #: A8424

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2797

Isotype: IgG

Swiss Prot: O43555

Purity: Affinity purification

For research use only.

Background:

This gene encodes a secreted peptide hormone and member of the gonadotropin-releasing hormone (GnRH) family of proteins. The encoded protein regulates reproductive function by stimulating the production and release of the gonadotropins follicle-stimulating hormone (FSH) and luteinizing hormone (LH). The encoded protein may inhibit endometrial, ovarian, prostate, and breast cancer cell proliferation. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).