

GTF3A

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:42kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human GTF3A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AP2; TFIIIA;

Catalog #:A8426

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2971

Isotype:IgG

Swiss Prot:Q92664

Purity:Affinity purification

For research use only.

Background:

The product of this gene is a zinc finger protein with nine Cys

To place an order, please [Click HERE](#).