

HLA-E

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 40kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human HLA-E

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MHC; QA1; EA1.2; EA2.1; HLA-6.2;

Catalog #: A8429

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3133

Isotype: IgG

Swiss Prot: P13747

Purity: Affinity purification

For research use only.

Background:

HLA-E belongs to the HLA class I heavy chain paralogues. This class I molecule is a heterodimer consisting of a heavy chain and a light chain (beta-2 microglobulin). The heavy chain is anchored in the membrane. HLA-E binds a restricted subset of peptides derived from the leader peptides of other class I molecules. The heavy chain is approximately 45 kDa and its gene contains 8 exons. Exon one encodes the leader peptide, exons 2 and 3 encode the alpha1 and alpha2 domains, which both bind the peptide, exon 4 encodes the alpha3 domain, exon 5 encodes the transmembrane region, and exons 6 and 7 encode the cytoplasmic tail.

To place an order, please [Click HERE](#).