

MAT2A

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 44kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MAT2A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MATA2; MATII; SAMS2;

Catalog #: A8436

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4144

Isotype: IgG

Swiss Prot: P31153

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene catalyzes the production of S-adenosylmethionine (AdoMet) from methionine and ATP. AdoMet is the key methyl donor in cellular processes.

To place an order, please [Click HERE](#).