

MEIS2

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 52kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MEIS2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MRG1; HsT18361;

Catalog #: A8437

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4212

Isotype: IgG

Swiss Prot: O14770

Purity: Affinity purification

For research use only.

Background:

This gene encodes a homeobox protein belonging to the TALE ('three amino acid loop extension') family of homeodomain-containing proteins. TALE homeobox proteins are highly conserved transcription regulators, and several members have been shown to be essential contributors to developmental programs. Multiple transcript variants encoding distinct isoforms have been described for this gene.

To place an order, please [Click HERE](#).