

MYL1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 21kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MYL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MLC1F; MLC3F;

Catalog #: A8438

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4632

Isotype: IgG

Swiss Prot: P05976

Purity: Affinity purification

For research use only.

Background:

Myosin is a hexameric ATPase cellular motor protein. It is composed of two heavy chains, two nonphosphorylatable alkali light chains, and two phosphorylatable regulatory light chains. This gene encodes a myosin alkali light chain expressed in fast skeletal muscle. Two transcript variants have been identified for this gene.

To place an order, please [Click HERE](#).