

NDUFA7

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 13kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NDUFA7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

B14.5a;

Catalog #: A8441

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4701

Isotype: IgG

Swiss Prot: O95182

Purity: Affinity purification

For research use only.

Background:

This gene encodes a subunit of NADH:ubiquinone oxidoreductase (complex I), which is a multiprotein complex located in the inner mitochondrial membrane. Complex I functions in the transfer of electrons from NADH to the respiratory chain.

To place an order, please [Click HERE](#).