

NRGN

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 8kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NRGN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RC3; hng;

Catalog #: A8444

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4900

Isotype: IgG

Swiss Prot: Q92686

Purity: Affinity purification

For research use only.

Background:

Neurogranin (NRGN) is the human homolog of the neuron-specific rat RC3/neurogranin gene. This gene encodes a postsynaptic protein kinase substrate that binds calmodulin in the absence of calcium. The NRGN gene contains four exons and three introns. The exons 1 and 2 encode the protein and exons 3 and 4 contain untranslated sequences. It is suggested that the NRGN is a direct target for thyroid hormone in human brain, and that control of expression of this gene could underlay many of the consequences of hypothyroidism on mental states during development as well as in adult subjects.

To place an order, please [Click HERE](#).