

POLG

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 140kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human POLG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PEO; MDP1; SCAE; MIRAS; POLG1; POLGA; SANDO; MTDP4A; MTDP4B;

Catalog #: A8451

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5428

Isotype: IgG

Swiss Prot: P54098

Purity: Affinity purification

For research use only.

Background:

Mitochondrial DNA polymerase is heterotrimeric, consisting of a homodimer of accessory subunits plus a catalytic subunit. The protein encoded by this gene is the catalytic subunit of mitochondrial DNA polymerase. The encoded protein contains a polyglutamine tract near its N-terminus that may be polymorphic. Defects in this gene are a cause of progressive external ophthalmoplegia with mitochondrial DNA deletions 1 (PEOA1), sensory ataxic neuropathy dysarthria and ophthalmoparesis (SANDO), Alpers-Huttenlocher syndrome (AHS), and mitochondrial neurogastrointestinal encephalopathy syndrome (MNGIE). Two transcript variants encoding the same protein have been found for this gene.

To place an order, please [Click HERE](#).