

SLCO1A2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 74kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SLCO1A2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OATP; OATP-A; OATP1A2; SLC21A3;

Catalog #: A8452

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6579

Isotype: IgG

Swiss Prot: P46721

Purity: Affinity purification

For research use only.

Background:

This gene encodes a sodium-independent transporter which mediates cellular uptake of organic ions in the liver. Its substrates include bile acids, bromosulphophthalein, and some steroidal compounds. The protein is a member of the SLC21A family of solute carriers. Alternative splicing of this gene results in multiple transcript variants.

To place an order, please [Click HERE](#).