

TBXA2R

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 37kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TBXA2R

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TXA2-R; BDPLT13;

Catalog #: A8458

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6915

Isotype: IgG

Swiss Prot: P21731

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the G protein-coupled receptor family. The protein interacts with thromboxane A2 to induce platelet aggregation and regulate hemostasis. A mutation in this gene results in a bleeding disorder. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).