

TGFB111

Reactivity:Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:50kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TGFB111

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HIC5; ARA55; HIC-5; TSC-5;

Catalog #:A8459

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7041

Isotype:IgG

Swiss Prot:O43294

Purity:Affinity purification

For research use only.

Background:

This gene encodes a coactivator of the androgen receptor, a transcription factor which is activated by androgen and has a key role in male sexual differentiation. The encoded protein is thought to regulate androgen receptor activity and may have a role to play in the treatment of prostate cancer. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).