

TGFB3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 47kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TGFB3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ARVD; RNHF; ARVD1; TGF-beta3;

Catalog #: A8460

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7043

Isotype: IgG

Swiss Prot: P10600

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the TGF-beta family of proteins. The encoded protein is secreted and is involved in embryogenesis and cell differentiation. Defects in this gene are a cause of familial arrhythmogenic right ventricular dysplasia 1.

To place an order, please [Click HERE](#).