

ZBTB17

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 88kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ZBTB17

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MIZ-1; ZNF60; ZNF151; pHZ-67;

Catalog #: A8463

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7709

Isotype: IgG

Swiss Prot: Q13105

Purity: Affinity purification

For research use only.

Background:

This gene encodes a zinc finger protein involved in the regulation of c-myc. The symbol MIZ1 has also been associated with PIAS2 which is a different gene located on chromosome 18.

To place an order, please [Click HERE](#).