

TAF15

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 62kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TAF15

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Npl3; RBP56; TAF2N; TAFII68;

Catalog #: A8465

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8148

Isotype: IgG

Swiss Prot: Q92804

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the TET family of RNA-binding proteins. The encoded protein plays a role in RNA polymerase II gene transcription as a component of a distinct subset of multi-subunit transcription initiation factor TFIID complexes. Translocations involving this gene play a role in acute leukemia and extraskeletal myxoid chondrosarcoma, and mutations in this gene may play a role in amyotrophic lateral sclerosis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

To place an order, please [Click HERE](#).