

GBF1

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 206kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GBF1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ARF1GEF;

Catalog #: A8468

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8729

Isotype: IgG

Swiss Prot: Q92538

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the Sec7 domain family. The encoded protein is a guanine nucleotide exchange factor that regulates the recruitment of proteins to membranes by mediating GDP to GTP exchange. The encoded protein is localized to the Golgi apparatus and plays a role in vesicular trafficking by activating ADP ribosylation factor 1. The encoded protein has also been identified as an important host factor for viral replication. Multiple transcript variants have been observed for this gene.

To place an order, please [Click HERE](#).