

F2RL3

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:41kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human F2RL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PAR4;

Catalog #:A8471

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9002

Isotype:IgG

Swiss Prot:Q96RI0

Purity:Affinity purification

For research use only.

Background:

Coagulation factor II (thrombin) receptor-like 3 (F2RL3) is a member of the large family of 7-transmembrane-region receptors that couple to guanosine-nucleotide-binding proteins. F2RL3 is also a member of the protease-activated receptor family. F2RL3 is activated by proteolytic cleavage of its extracellular amino terminus. The new amino terminus functions as a tethered ligand and activates the receptor. F2RL3 is activated by thrombin and trypsin.

To place an order, please [Click HERE](#).