

INADL

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 196kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human INADL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Cipp; PATJ; hINADL; InaD-like;

Catalog #: A8476

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10207

Isotype: IgG

Swiss Prot: Q8NI35

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein with multiple PDZ domains. PDZ domains mediate protein-protein interactions, and proteins with multiple PDZ domains often organize multimeric complexes at the plasma membrane. This protein localizes to tight junctions and to the apical membrane of epithelial cells. A similar protein in Drosophila is a scaffolding protein which tethers several members of a multimeric signaling complex in photoreceptors.

To place an order, please [Click HERE](#).